
R2U-TOX-ASSAY Human Diversity in a 96 well Format for Improved Drug Discovery





On-demand hiPSC-based cell models in *ready-to-use* 96 well plate format for phenotypic and toxicological screening

The combination of human induced pluripotent stem cell (hiPSC) technologies with our advanced vitrification-based cryotechnology platform leads to:

- Time and resources savings due to a more reliable design of screening campaigns
- Stock keeping of mature hiPSC-based model systems
- High-quantity production and cryopreservation of adherent cell models using scalable platform technologies
- Reduction of animal testing by using human model systems with higher predictive value
- Facilitating early-killer experiments improving efficiency in drug discovery

Ultra-low-temperature stable 96 well plate enables ice-free cryopreservation (vitrification) of adherent cell models






Neural R2U plate

Neural co-culture consisting of astrocytes & midbrain neurons

Non-frozen control




Vitrified product post warming




Preservation of neural network morphology and cell number

Transmitted light images of a vitrified (24h after warming) and non-frozen neural co-culture.

Non-frozen control




Vitrified product post warming



Maintenance of neural-specific gene and protein expression levels (ICC, qPCR, FACS)

Key biomarkers: GFAP, MAP2, TUBB3, S100b

Fluorescence images of a vitrified (24h after warming) and non-frozen neural co-cultures visualized using the markers for astrocytes (GFAP) and neurons (MAP2).



Neural co-culture after Fluo-4 loading.



Vitrified product post warming



Positive action potential measured by relative fluorescence signal in comparison between a non-frozen and vitrified cell (24h after warming).





Fraunhofer Institute for Biomedical Engineering (IBMT) Dr. Julia Neubauer Joseph-von-Fraunhofer-Weg 1 66280 Sulzbach Germany

julia.neubauer@ibmt.fraunhofer.de www.ibmt.fraunhofer.de

